




SVC-BASED STRESS TEST APPLICATION (USER GUIDE)




USER INTERFACE





PROCEDURE (ANALYZE)


- *Please load your SVC Model before attempting to analyze your data
- 1. Complete all entries with your data
- 2. The data analysis will be based on your loaded SVC Model
- 3. Press the 'Run' button to start analyze your data (only available if SVC Model has been loaded)
- 4. The analysed result will be displayed as either 'Yes' or 'No'
- 5. Press the 'Reset' button to analyze your next data






PROCEDURE (MODEL - LOADING SVC MODEL)

- 1. Enter your SVC Model address by manually typing it into the entry or;
- 2. Press the 'Browse File' button to help you locate your SVC Model and autofill the entry address *If you have trained your SVC Model by using the app, your SVC Model should be located in the app folder under 'models' (Refer Pg. 6)
- 3. Press the 'Load' button to load your selected SVC Model
- 4. The current loaded SVC Model will be displayed on the 'Analyze' page









PROCEDURE (MODEL – TRAINING SVC MODEL)

- 1. Enter your .csv dataset file by manually typing it into the entry or;
- 2. Press the 'Browse File' button to help you locate your .csv dataset file and autofill the entry
- 3. Enter your preferred file name for the new SVC Model
- 4. Press the 'Train' button to start training the SVC Model
- 5. Your newly trained SVC Model will be located in the app folder under 'models'








DATASET FOR SVC MODEL TRAINING

- 1. The app only accepts datasets in .csv file format
- 2. If you want to create a new dataset file, ensure that the .csv data format is as shown:





APPLICATION UPDATES

- 1. Visit our GitHub page to download the latest app new updates or releases at https://github.com/durio-dev/svc-stress-app.git
- 2. Locate and delete all files in your current app folder
- 3. Unzip the downloaded 'svc-stress-app-main.zip' and copy the content of all files under 'releases' and paste them into your current app folder



THANK YOU